NSK


Multi Function Ultrasonic Scaler

Varios 370

Varios 370
iPiezo engine

OPERATION MANUAL

Please read this Operation Manual carefully before use,
and file for future reference.

Rx Only

OM-E0620E 001 **MADE IN JAPAN**

◆ Original Operation Manual**◆ Classifications of equipment**

- Type of protection against electric shock:
 - Class II equipment 
- Degree of protection against electric shock:
 - Type BF applied part: 
- Method of sterilization or disinfection recommended by the manufacture:
 - See 12. Sterilization
- Degree of protection against ingress of water as detailed in the current edition of IEC 60529:
 - Foot Control: IPX1 (Protected against vertically falling water drops)
- Degree of safety of application in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide:
 - EQUIPMENT not suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide.
- Mode of operation:
 - Continuous operation

Intended to Use

This product is intended only for dental clinic /dental office use. This device generates ultrasonic waves intended for use in dental applications such as scaling, root canal treatment, periodontal and cavity preparation.

Rx Only Caution: U.S. Federal law restricts this device to sale by or on the order of a dentist.

**Cautions for handling and operation**

- Read these cautions carefully and use only as intended or instructed.
- Safety instructions are intended to avoid potential hazards that could result in personal injury or damage to the device. Safety instructions are classified as follows in accordance with the seriousness of the risk.

**WARNING**

- A hazard that could result in bodily injury or damage to the device if the safety instructions are not followed.

**CAUTION**

- A hazard that could result in light or moderate bodily injury or damage to the device if the safety instructions are not followed.

NOTICE

- General information needed to operate the device safely.

**WARNING**

- TO PREVENT ELECTRIC SHOCK Do not unplug the AC Adaptor with wet hands.
- TO PREVENT ELECTRIC SHOCK Make sure not to get water on the Control Unit.
- TO PREVENT ELECTRIC SHOCK Do not touch the handpiece backend electrical connections.
- TO PREVENT ELECTRIC SHOCK Use an electrical outlet that is grounded.
- If you feel any abnormality such as vibration, heat generation, abnormal noise, etc., prior or during the use of the unit, stop using it immediately.
- This product is Medical Electrical equipment Electromagnetic compatible (EMC). As described in the accompanying documentation.
- Portable and mobile RF communications equipment can affect Electrical Medical equipment. Do not use RF equipment in close proximity to the product.
- When installing the product, provide space of approximately 10cm around the Control Unit for easy access to the inlet and the AC Adaptor.

- Be sure to always use with sufficient coolant water. Insufficient coolant water could cause a burn of the affected area or the overheating of the instrument.
- Be aware of the tip end position. Constantly move the tip so that it does not stagnate at one point. If the tip is stagnated at a single point, there will not be enough cooling, and this could result in overheating of the tooth or root canal, or excessive cutting of the tooth or root canal.
- USE ONLY NSK genuine Tips when using NSK Varios Ultrasonic Scaler (Varios 370 or Varios 370 Lux) problems such as damage, failure and accident of Handpieces resulting from use of Non-NSK Tips are not included in the warranty. The following are the possible failure that could happen when using the Non-NSK Tips;
 - Vibration failure caused by using non conforming screws.
 - Patients accidental ingestion of broken Tips.
 - Damage of thread ridge of handpiece.
- You must use the Tip within the power range described on the Tip Power Guide. If you use it out of the power range, the Tip might break or damage an operative site.
- When operating the product always consider the safety of the patient.
- Use by medical professional, such as doctor or dental hygienist, is intended.
- Check the vibration outside the patient's oral cavity before use. If any abnormalities are found, stop using immediately and contact dealer.
- Do not drop, hit, or give excessive shock on the Control Unit/Handpiece. This could lead to electric shock or product failure. In case of giving shock to an optic handpiece, illumination parts that are equipped with the handpiece could be damaged and drop out of the handpiece.
- To prevent possible tooth plane damage and handpiece overheating, always use with sufficient water.
- Do not sterilize by ultraviolet light. Handpiece could discolor.
- Sterilize the Tip (except for all diamond coated tip), Handpiece, Tip Holder, Tip Cover S and Tip Wrench by autoclaving. Wipe the Control Unit, Tip Holder, Tip, AC Adaptor, Foot Control, Tip Card, Tip Power Guide, and Handpiece Cord including the cover with alcohol soaked cotton.
- If chemical, solvent or antiseptic solution is deposited on this product, immediately wipe it away. Discoloration or deformation may occur if left.
- Do not disassemble or alter the handpiece/Control Unit.
- Keep away from patients with cardiac pacemakers.
- Keep away from explosive substances and flammable materials. Do not use for patients anesthetized under laughter gas. (Nitrous oxide)
- Staring directly into LED may result in eye damage. Do not turn it to the eyes of the patient or operator.
- This product needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information.
- Electromagnetic Compatibility: The use of ACCESSORIES, transducers and cables other than those specified, with the exception of transducers and cables sold by the manufacturer of this product as replacement parts for internal components, may result in increased EMISSIONS or decreased IMMUNITY of this product.
- This product should not be used adjacent to or stacked with other equipment and that if adjacent or stacked use is necessary, this product should be observed to verify normal operation in the configuration in which it will be used.
- If any water drops remain on the handpiece or handpiece cord after autoclaving, wipe them off. Staining may result if left.
- Users are responsible for operation control and maintenance.
- Do not use water tubes of other companies. Doing so may cause a malfunction.
- There is the judgment that applies this product to a patient in the user side.
- Do not apply excessive power to the Tip. It may damage the teeth because of the ultrasonic vibration.



CAUTION

- During operation, high frequency oscillations in the handpiece and handpiece cord may affect computer and LAN Noise may be heard during operation near a radio receiver.
- Be sure to turn off the Power/Volume Knob after use. Remove the AC Adaptor when not using it for long periods of time.
- The surface of the Control Unit may feel warm after prolonged use. Turn the Power/Volume Knob to OFF and cool the device.
- Users are responsible for operational control, maintenance and inspection.
- Clean/sterilize the product immediately after using it. Then store it. Leaving it non-sterile might lead to failure.
- When you have not used the product for long time and use it again, check the operation before use.
- When abnormalities are found with a Control Unit and /or an AC adaptor, pull the AC adaptor from AC Outlet immediately.
- Do not turn the Power/Volume Knob to ON or OFF needlessly.
- No special training is required for this device.
- Applied parts for patient and/or operator are/is Tip and Handpiece.
- When using a file as a tip for root canal treatment (for enlarging, cleaning, etc.), use it after confirming its characteristics and areas where it can be applied.
- Operators and all others in the area must wear eye protection and a mask when operating this product.
- Before use, be sure to always check that the power level is set to the proper settings.
- During vibration, do not grasp the tip nor make contact with the tip for an extended time.
- Do not use tips that are damaged, bent, corroded, or show other imperfections. The tip could break or fail to vibrate.
- Do not apply more pressure than required to the tip end. This could cause the tip to break or warp. Do not twist or strain the tip, and do not use it as a replacement for surgical tools. This could cause the tip to deform or break.
- If the tip cracks or breaks during use, immediately stop operation, and replace with a new tip. Operation with a damaged tip could injure the gums, periodontal tissue, or other locations.
- The diamond coating on the tip surface may peel off depending on the tooth condition or root canal condition. If the diamond coating peels off, the cutting efficiency becomes less, resulting in reduced cutting performance. If you sense this type of condition, replace with a new tip.
- All sterilizable components of this device are delivered in a non-sterile condition and must be autoclaved prior to use. Refer to "4. Sterilization before first time use".
- The diamond coated tips are "Single use". Do not reuse.
- U.S. Federal law restricts this device to sale by or on the order of a licensed physician.

* Operation Principle

A sinusoidal electrical signal, at ultrasonic frequency ($f > 20\text{kHz}$), is delivered by a generator. This signal is applied to the 'piezoelectric ceramic' located inside the transducer. Piezoelectric ceramic converts this signal into mechanical vibrations. These vibrations are at the same ultrasonic frequency as the electrical signal. The mechanical vibrations are propagated towards the distal end of the transducer. The "TIP" insert, which is attached at the distal end of the transducer, vibrates at ultrasonic frequencies and makes it possible to achieve the aimed purpose.

1. Component Names

No.	Parts Name	Quantity
1	Control Unit (with Handpiece Cord Unshielded 2M)	1
2	Handpiece (Varios2 or Varios2 Lux)	1
3	Control Unit Holder	1
4 ^{*1}	AC Adaptor (Unshielded cord 5M)	1
5	Handpiece Holder	1
6	Foot Control (Unshielded 4M)	1
7	Water Tube Set	1
8	Tip (G4,G6, G8)	1

No.	Parts Name	Quantity
9	Rubber Pad	4
10	Tip Wrench	1
11	Spanner Wrench (5x8)	2
12	O Ring	2
13	Tip Cover S (Option)	-
14	Tip Holder (Option)	-
15	Set Screw ^{*3}	3


*¹ Adaptor will change according to region

*² Only 120 V


*³ These are not on photo above.

2. Name and Function of each part

Front side (Control Unit With Control Unit Holder)

Back side

3. Prior to Operating System


3-1 Set the AC Adaptor

Insert each plug into appropriate connector.

① Insert the Plug Head into the Adapter Body as shown in the right figure.

② Slide the Plug Head on to the Adapter Body.

To release, push the Release Button shown on the right figure and remove the Plug Head from the Adapter Body.

3-2 Connecting

Insert each plug into appropriate connector. (Fig.1)

- 1 Align the "△" Mark on the Foot Control Connector and Foot Control Plug and connect those firmly into Foot Control Connector.
- 2 Connect Water Tube (Non-Filter Side, refer to Fig. 23 for detail) firmly into Water Tube Connector.
- 3 Connect AC Adaptor into DC Connector.


Fig.1

CAUTION

- Insert each plug and tube firmly into the connector. If connection is incomplete, it may cause a malfunction.
- Ensure Power is OFF on the Control Unit during the AC Adaptor Connection.
- Do not connect the cord in wall outlet before connecting DC Connector.
- Do not pull the AC Adaptor forcibly.
- Do not disconnect the AC Adaptor while pressing on the Foot Control.
- Turn OFF the power to connect or disconnect the cords and plugs.

3-3 Disconnecting

3-3-1 DC Plug and Foot Control Plug

Simply pull out plugs from the Control Unit.

3-3-2 Water Tube (Fig.2)

Depress the white ring to remove the water tube.


Fig.2

CAUTION

- Shut off the main valve of the water injection, and then remove the tube.
- It requires the water removal before the water tube disconnection.

3-4 Handpiece Holder

The handpiece holder attached to the Control Unit keeps the handpiece in place when it is not in use. To mount the Handpiece Holder on a vertical surface follow the procedures below.

- 1) Make a 3 mm hole in the Double-Face tape behind the screw hole in handpiece holder. (Fig. 3)
- 2) On the vertical surface make a 3.2~3.5 mm hole in the plate where the holder will be mounted. Align the hole with the screw hole on the holder. Make certain the surface is clean then adhere the holder to the plate, and attach. Insert the Set Screw through the two holes on the holder and the plate. Fasten the Set Screw with the Nut from the rear of the plate. (Fig. 4)


Fig. 3


Fig. 4

3-5 Handpiece Cord Holder and Control Unit Holder

The Control Unit can be attached to the chair unit using the control unit holder.

- 1) Make two 3.2~3.5mm holes in the mounting surface. Insert the set screw in the hole from the holder side and secure it with the nut.
- 2) Align the chaser and Slide the Control Unit into the Control Unit Holder.
- 3) Tuck the Handpiece Cord into Handpiece Cord Holder.


Fig. 5

- Do not pull or bend the handpiece cord forcibly as it might cause damage to the handpiece, and result in improper water flow.
- You can mount Control Unit Holder on a top surface or bottom surface.

3-6 Rubber Pad

To prevent the control unit from slipping apply rubber pads.

- 1) Clean bottom of the Control Unit.
- 2) Fit the Rubber Pad appropriate place as shown on the Fig. 6.


Fig. 6

CAUTION

- Control Unit Holder is not necessary when rubber pads are applied.

4. Sterilization before first time use

NSK Varios components are not sold sterile and should be sterilized prior to first time you use. Autoclave sterilization is recommended.

Please note the following when sterilize components.

Follow your local directives, standards, and guidelines.

Sterilize components in a U.S. FDA-cleared sterilization pouch, complaint to ISO 11607-1.

Ensure that you only remove dry sterile goods from the sterilizer.

Store sterile goods dust-free and dry.

NOTE

- NSK recommends sterilization according to ANSI/AAMI ST79 and EN 13060 (class B).
- Follow the manufacturer's instructions for use that came with your sterilizer.
- For more information about sterilization standards, refer the medical device regulations specific to your location.

■ Autoclave Procedure

- 1) Pick a tip up from the tip case. Take handpiece out of the packing bag before sterilization.
- 2) Insert sterilizable components into a U.S. FDA-cleared sterilization pouch, complaint to ISO 11607-1. Seal the pouch.
- 3) Autoclave conditions are below:

< When sterilizing a Tip >

Type	Steam Sterilization Cycle					
	Gravity Displacement		Pre-Vacuum (Dynamic Air Removal)			
Minimum Temperature	132°C	135°C	132°C	134°C*	135°C	
Full Cycle Time	15 minutes	10 minutes	4 minutes	3 minutes	3 minutes	
Minimum Drying Time	15 minutes	30 minutes	20 minutes	30 minutes	16 minutes	

* Not for use in the U.S.

<When sterilizing components except for a Tip >

Type	Steam Sterilization Cycle			
	Gravity Displacement			Pre-vacuum
Minimum Temperature	121°C	132°C	135°C	135°C
Full Cycle Time	20 minutes	15 minutes	10 minutes	10 minutes
Minimum Dry Time	10 minutes	10 minutes	10 minutes	15 minutes

- Do not exceed 138°C during sterilization.

CAUTION

- 4) Keep the components in the autoclave pouch until ready to use.
- 5) Store in a dry and dust-free location up to the period specified by the sterilization pouch manufacturer. If sterility cannot be confirmed, sterilize again prior to use.

5. Mounting and Removing the Handpiece

Align the Dots on the Handpiece and the Handpiece Cord. Push handpiece into connector.

To detach the handpiece, grip the Handpiece and plug of the handpiece cord and separate in a straight motion. (Fig. 7)


Fig. 7

- To avoid Electrical Shock Do not touch the handpiece backend electrical contacts.

WARNING

- Always confirm that the handpiece is correctly seated and locked into place.
- Do not use a Handpiece other than that supplied with this device.

CAUTION**6. Mounting and Removing Tip**

- 1 Turn TIP lightly by hand, and install it.
- 2 Tip will insert from the bottom hole of Tip Wrench. Align the four corner of the Tip base area into the four corner of Tip Wrench. And turn it clockwise until it clicks.
 - * Do not touch the top part of TIP to avoid an injury.
- 3 To remove the Tip, turn it counterclockwise with the Tip Wrench.



Fig. 8

**CAUTION
for Tip
Usage**

- Check the Tip before use. (Flush, Damage, Bending or Rust)
- Do not exceed Maximum Power Level for Tip. Damage to tooth structure and Tip may result.
- Do not hit metal or prosthetic crown etc. except when removing them. Tip could break and fall into mouth.
- Do not hit gingival, mucosa and/or skin directly. It could cause damage and/or burn injury.
- During its operation, the whole tip vibrates due to ultrasonic waves. Do not bring any part of the tip into contact with soft tissue, gums or skin. Regardless of use of coolant water, the tip in vibration could cause burn injury.
- The “dry tip* ” (a tip which does not require coolant water) utilizes heating to some extent for effective treatment. During or after its operation, the whole of such tip could be overheating. Do not bring any part of the tip into contact with soft tissue, gums or skin. Ensure that the power setting of the Control Unit does not exceed the maximum power written on the Tip Case label. Adjust the power from the lowest level. During use of the tip, ensure that burn injury does not occur in or around the treatment area. Also ensure that the tip does not affect the treatment area and oral cavity.
* E5, E6, E7, E8, G21, G22, G28
- Do not sharpen and/or bend the Tip. Tip may damage and not generate enough vibration during scaling.
- During cutting, Tip will gradually wear away, as the Tip wears the stroke will get smaller and decrease cutting efficiency. When level drops too far, change the Tip. (Tip Card check)
- DO ENSURE when securing Tip to use the Tip Wrench as supplied, inefficient cutting will result.
- DO ENSURE before attaching Tip, Cleanliness of the tip screw, inefficient cutting will result.
- To avoid personal injury ENSURE Tip is removed prior to disconnecting the handpiece from the handpiece cord.
- If you feel the Tip is not vibrating, remove it from an operative site, and press the Foot Control again. If this does not improve the condition, Ensure the Tip is secure, turn the power off and restart it.
- When mounting the Tip, always use grooves and Tip Wrench as supplied.
- Turn water off when using Tip that does not require water.
- Tip Wrench is consumable For reliable operation replace annually.

7. Operating Procedures

7-1 Power On (Fig. 9)

Connect the AC Cord to the wall outlet. Rotate the Power/Volume Knob on the Control Unit. (Power indicator will light.)


Fig. 9

7-2 Power Level Setting

DO ENSURE Power setting does not exceed the recommended Power Level (Tip-Power Guide or Tip case label (Fig.10) included in the package.)

Set the power level with the Power/Volume Knob on the Front Panel. Make sure the power level is set in the appropriate range for the attached Tip.

* Power Level for each mode



Fig. 10

- NOTICE**
- Rotate the Power/Volume Knob to increase or decrease the Power Level.
 - If the Power Level is set to 0 (zero), Tip will not oscillate but water comes out from the handpiece.

7-3 Operate Varios 370 / 370 Lux

Tip vibration will begin when the Foot Control is depressed. Also, Output Indicator will light. (For Varios2 Lux, Handpiece LED will illuminate.)

7-3-1 Water Supply Volume Adjustment

Turn the Water Adjustment Knob clockwise gradually to increase the supply volume. (Fig. 11)

◆ During the Handpiece operation :

Possible: Power Level and Water Volume adjustment.


Fig. 11

7-4 After the Treatment

Release the Foot Control and Power off the Control Unit. Close the dental unit's water valve.

◆ How to drain the inside of the Handpiece

- 1) Shut off the water main valve.
- 2) With the Tip and Handpiece attached, turn the Power/Volume Knob to its lowest setting and then continue stepping on the Foot Control until water stops coming from the Tip.
- 3) Remove the Tip and Handpiece.

NOTICE • LED of the handpiece will remain 'On' for approx 5 seconds after Foot Control is released. (Varios2 Lux)

7-5 Protection Circuit

The Varios 370 has an internal Protection Circuit to prevent overheating. If the internal temperature of the Varios 370 becomes hot, such as when the device is used for more than 10 minutes with the power level set higher than G7, the output of the device will be reduced automatically to G7. The device will continue to operate at this reduced level until the internal temperature has reduced to below the protection circuit threshold. During this period, the maximum output setting that can be used is G8 and to increase output from G7 to G8, output must first be reduced to G5 or lower, then increased to the maximum of G8. Maximum output can be restored by turning the device off for 15 minutes or longer.

NOTICE • During Protection Circuit function, the Control Unit can not increase the Power Level more than 8.

8. Provided Scaler Tips

G4

The end of the Tip is thin and for supragingival fine scaling and interdental scaling. The round cross-section allows tooth surfaces to be finished without causing damage.

Apply the top of the Tip on the tooth plane and move it sideways finely in the same way as G8 Tip. (Fig. 12)

Fig. 12

G6

Removal of supra and subgingival calculus. It provides easy access to interdental spaces and narrow pockets.

Insert the top of the Tip into the periodontal pocket and move it slowly. The top of the Tip is sharp so that it could remove tartar on long coroner and retracted gingival. (Fig. 13)
Clean periodontal pocket at low power. (Set the level less than "Power 5" at P mode.)


Fig. 13

G8

Removal of supragingival and interdental calculus. This Tip can be used in all quadrants and is very useful for the removal of hard calculus.

Apply the top of the Tip on the tooth plane and move it sideways finely along the neck of tooth. (Fig. 14)

Fig. 14


CAUTION

- Tip is article of consumption. We recommend periodical replacement. About time of replacement, check the Tip Card.
- Additional Tips and accessories are available for order in the catalog.


◆ How to use the Tip Card

- 1) Place the neck of the Tip in the cut out.
- 2) Check wear of the Tip.
- 3) Reference the green, yellow and red line to check wear of the Tip. *Ref. Fig. 15. At NSK we recommend to replace a Tip when the Tip meets the yellow line (wear of 1mm) to guarantee safe and effective use.

CAUTION

- Tips are consumables. The efficiency of dental scaling decreases approximately 25% when the top of the Tip wears 1 mm and approximately 50% when it wears 2 mm. In addition, the vibration condition changes according to the wear, which may damage a patient's tooth surface. Check the Tip wear condition with the Tip Card periodically, and replace the Tip with a new one in good time.

9. How to Use Tip Cover S (Option)

Grip the Tip Cover S and insert it to the Tip.
To remove, grip the Tip Cover S and the handpiece & pull.
(Fig. 17)


Fig. 17

CAUTION

- Carefully insert the Tip into the Tip Cover S. Avoid injuring the fingers.
- The Tip Cover S is not designed for use as a Tip changing tool.

10. Holder

10-1 Handpiece Holder

While the Handpiece is not in use, put the Handpiece on the Handpiece Holder. (Fig. 18)

NOTICE • To prevent injury, always mount Scaler Tip Cover S.


Fig. 18

10-2 Tip Holder (Option)

The Tip Holder is Autoclavable and hold up to 5 tips at once. To Autoclave, tilt the tips in the direction of the arrow in Fig. 19.


Fig. 19

11. Care and Maintenance

11-1 Cleaning of Optic Fiber (Varios2 Lux)

Wipe the debris off the end of the Optic Fibers at the handpiece with alcohol. (Fig. 20)

NOTICE • Do not use any sharp pointed tools to clean the Optic Fiber End Face. In case the light degradation, contact dealer.


Fig. 20

11-2 Changing O-Ring

Handpiece Cord

An O-Ring is located in the Handpiece Cord Connector. Replace O-Ring in the groove. (Fig. 21)

* Optional O-Ring: Order Code D0310020080

Fig. 21

11-3 Changing Water Filter (Option)

Change the Water Filter as it may necessary.

- 1) Close the water valve of the Chair Unit.
- 2) Mount two Spanner Wrenches (5x8) and turn those as shown in Fig. 22.
- 3) When the Water Filter case is separated, the Water Filter can be removed as shown in Fig. 23.
- 4) Replace with new (Order Code U387042) and reassemble the filter in the reverse order.


Fig. 22


Fig. 23

12. Periodical Maintenance Checks

In order to use the product properly and safely, perform operation checks once a year following the steps shown in "3. Prior to Operating System", "5. Mounting and Removing the Handpiece", "6. Mounting and Removing Tip", and "7. Operating Procedures". Contact your Authorized NSK Dealer if any abnormalities are found.

13. Sterilization

- All sterilizable components of this device must be sterilized after each patient.
- The recommended method of sterilization is by autoclave sterilization.
- ONLY the Tip (except diamond coated scaler tip), Handpiece, Tip Wrench, Tip Holder and Tip Cover S can be autoclaved. The diamond coated Scaler Tip is to be sterilized prior to first use, then disposed of as this Tip is single use only.

- NOTE**
- NSK recommends sterilization according to ANSI/AAMI ST79 and EN 13060 (class B).
 - Follow the manufacturer's instructions for use that came with your sterilizer.
 - For more information about sterilization standards, refer the medical device regulations specific to your location.

The surface of the handpiece can be washed and disinfected with a Thermo-Disinfector prior to Sterilization.

■ Autoclave Procedure

- 1) Remove the Tip after use. (Refer to 6. Mounting and Removing Tip)
- 2) Clean all the surfaces of the components using a small tooth brush with isopropyl alcohol to remove dirt and debris then dry thoroughly.
- 3) Insert sterilizable components into a U.S. FDA-cleared sterilization pouch, complaint to ISO 11607-1. Seal the pouch.
- 4) Autoclave conditions are below:

< When sterilizing a Tip >

Type	Steam Sterilization Cycle				
	Gravity Displacement		Pre-Vacuum (Dynamic Air Removal)		
Minimum Temperature	132°C	135°C	132°C	134°C*	135°C
Full Cycle Time	15 minutes	10 minutes	4 minutes	3 minutes	3 minutes
Minimum Drying Time	15 minutes	30 minutes	20 minutes	30 minutes	16 minutes

* Not for use in the U.S.

<When sterilizing components except for a Tip>

Steam Sterilization Cycle				
Type	Gravity Displacement			Pre-vacuum
Minimum Temperature	121°C	132°C	135°C	135°C
Full Cycle Time	20 minutes	15 minutes	10 minutes	10 minutes
Minimum Dry Time	10 minutes	10 minutes	10 minutes	15 minutes

- Do not exceed 138°C during sterilization.

CAUTION

- 4) Keep the components in the autoclave pouch until ready to use.
- 5) Store in a dry and dust-free location up to the period specified by the sterilization pouch manufacturer. If sterility cannot be confirmed, sterilize again prior to use.

CAUTION

- Do not sterilize by ultraviolet ray. The handpiece could discolor.
- If autoclaved with other instruments stained with chemical solution, it could strip the plating and make the surface black.
- Do not autoclave any parts (the Control Unit, AC Adaptor, Foot Control, Handpiece Cord, Tip Card, Tip Power Guide, O-Ring). Other than those that can be subjected to autoclave sterilization. Perform alcohol disinfection to the Control Unit, AC Adaptor, Foot Control, Handpiece Cord, Tip Card, Tip-Power Guide including after every patient.
- Do not wipe with, or clean or immerse in, high acid water or sterilizing solutions.
- Diamond coated Scaler Tip is for single use only. Do not sterilize and use again.
- Keep the product in suitable atmospheric pressure, temperature, humidity, ventilation, and sunlight. The air should be free from dust, salt and sulphur.
- Do not touch the product immediately after autoclaving may cause injury.

14. Troubleshooting

When trouble is found, please check the followings prior to consulting your dealer.

Problem	Probable Cause	Cause	Solution
No / Poor vibration	With power on and no Power Indicator	The AC Adaptor or the DC Plug is disconnected	Correctly insert the AC Adaptor or the DC Plug
	The Tip does not generate vibration, in spite of depressing the Foot Control	The Tip is not tightened firmly	Tighten the Tip until the Tip Wrench clicks
		Worn Tip	Replace the Tip
		Power has not been correctly adjusted for the Tip	According to Power Guide or Tip case label. Do not exceed recommendations.
		The Foot Control is disconnected	Connect the Foot Control correctly
		Failure of vibrator in the handpiece	Contact dealer*
Failure of internal components of the Foot Control	Contact dealer*		
The Tip is bent or broken	—	Power has not been properly adjusted for the Tip	According to Power Guide or Tip case label. Do not exceed recommendations.
The Tip is flying away	—	The Tip is not tightened firmly	Tighten the Tip until the Tip Wrench clicks
Noise from the handpiece	—	Power has not been properly adjusted for the Tip	According to Power Guide or Tip case label. Do not exceed recommendations.
		The Tip is not tightened firmly	Tighten the Tip until the Tip Wrench clicks
		Failure of vibration in the handpiece or the Control Unit	Contact dealer*
The handpiece is overheating	—	Power has not been properly adjusted for the Tip	According to Power Guide or Tip case label. Do not exceed recommendations.
		The Tip is not tightened firmly	Tighten the Tip until the Tip Wrench clicks
		Failure of vibration in the handpiece or the Control Unit	Contact dealer*
No / Poor water	The water does not reach to the Control Unit	—	Check the water circuitry and supply to the Control Unit. Water pressure : 0.1-0.5MPa (1-5kgf/cm ²)
	Check to see if water reaches the Control Unit	The Water Adjustment Knob is closed.	Turn the Water Volume Knob and adjust to the appropriate volume
		Disconnected water supply at low volume range. (less than 10ml/min.)	Adjust the Water Volume Knob and increase the Irrigation volume
		The Water Filter is clogged	Replace with new Water Filter (Refer to 11-3 Changing Water Filter (Option))
Water leakage	Water is leaking from the joint between the handpiece and the cord	O-Ring at the handpiece cord is worn or damaged	Replace with new O-Ring (Refer to 11-2 Changing O-Ring Handpiece Cord)
Attachment of the Control Unit Holder is loose	—	Holder worn	Contact dealer*

Problem	Probable Cause	Cause	Solution
Handpiece LED does not illuminate. (Varios2 Lux)	Tip oscillates, but Handpiece LED turns on and off	The handpiece is not connected into the Handpiece Cord correctly	Firmly insert the handpiece into the Handpiece Cord completely
Loss of the power output without operation	Power output is set 8 at G	Protection Circuit is activated	Power output will weaken automatically while continuous operation is over 10min at the setting of Maximum power at G mode. Releasing the foot from the Foot Control. Decrease the Power less than 5, once then increase the power again. (Refer to 7-5 Protection Circuit)

* Repairs cannot be made by the customer.

15. Specification

Type	VA370 Lux/VA370
Power source	AC 100 - 240 V 47 - 63Hz
Vibration Frequency	28 - 32 kHz
Maximum Output	11 W
Rated Power	25VA
Lighting	Varios 370 : No Varios 370 Lux: Yes
Dimension	W 80 x D 115 x H 32 mm (Body without knob and cord)
Weight	0.43 kg (Except Attachment)
Use Environment	Temperature 0 to 40 °C (The liquid must not freeze up) Humidity 30 to 75 % Atmospheric Pressure 700 to 1060 hPa
Store Environment	Temperature -10 to 60 °C Humidity 10 to 85 % Atmospheric Pressure 500 to 1060 hPa

16. Spare Parts

Model	Products	Order code	Model	Products	Order code
Water Tube Set		U1007080	Tip Holder	 	Z221A080
Water Supply Connector		U387030	Tip Cover S	 	Z217851
Water Filter		U387042	O-Ring		D0310020080
Spanner Wrench (5x8) X 2 pcs		Y1001301	Set screw X 3 pcs Rubber Pad		U1007091
Tip Wrench (CR-10)	 	Z221076	Double-Face Tape (For Handpiece Holder)		20002544

 Autoclavable up to max 135°C.

17. Disposing product

Consult with dealer from whom you purchased it about waste disposal.

18. Warranty

Manufacturer warrants its products to the original purchaser against defects in material and workmanship under normal practices of installation, use and servicing. Such expendable items as O-Rings are not covered by this warranty.

Symbols

TUV Rhineland of North America is a Nationally Recognized Testing Laboratory (NRTL) in the United States and is accredited by the Standards Council of Canada to certify electro-medical products with Canadian National Standards.

Follow the waste of electric and electronic equipment (WEEE) Directive (2012/19/EC) to dispose of the product and accessories.

Consult operation instructions.

Manufacturer.

Class II Equipment.

This conforms to CE European Directive of "Medical equipment directive 93/42/EEC."

Type BF applied part.

Authorised representative in the European community.

Protected against vertically falling water drops.

Autoclavable up to Max. 135°C.
*for detail see Sterilization.

The surface of the product can be washed and disinfected with a Thermo-Disinfectant prior to Sterilization.

Marking on the outside of Equipment or Equipment parts that include RF transmitters or that apply RF electromagnetic energy for diagnosis or treatment.

Single use. Do not reuse.

Rx Only Rx Only Caution: U.S. Federal law restricts this device to sale by or on the order of a dentist.

Guidance and manufacturer's declaration - electromagnetic emissions		
The Varios 370 / Varios 370 Lux is intended for use in the electromagnetic environment specified below. The customer or the user of the Varios 370 / Varios 370 Lux should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR11/EN55011	Group 1	The Varios 370 / Varios 370 Lux uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR11/EN55011	class B	The Varios 370 / Varios 370 Lux is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supply network that supplies buildings used for domestic purposes.
Harmonic emissions EN/IEC61000-3-2	class A	
Voltage fluctuations/flicker emissions EN/IEC61000-3-3	Complies	

Guidance and manufacturer's declaration - electromagnetic immunity			
The Varios 370 / Varios 370 Lux is intended for use in the electromagnetic environment specified below. The customer or the user of the Varios 370 / Varios 370 Lux should assure that it is used in such an environment.			
Immunity test	EN/IEC60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) EN/IEC61000-4-2	±6kV contact ±8kV air	±6kV contact ±8kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst EN/IEC61000-4-4	±2kV for power supply lines ±1kV for input/output	±2kV for power supply lines ±1kV for input/output	Mains power quality should be that of a typical commercial or hospital environment.
Surge EN/IEC61000-4-5	±1kV line(s) to line(s) ±2kV line(s) to earth	±1kV line(s) to line(s) ±2kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines EN/IEC61000-4-11	<5% Ut (>95% dip in Ut) for 0.5 cycle 40% Ut (60% dip in Ut) for 5 cycles 70% Ut (30% dip in Ut) for 25 cycles <5% Ut (>95% dip in Ut) for 5 secs	<5% Ut (>95% dip in Ut) for 0.5 cycle 40% Ut (60% dip in Ut) for 5 cycles 70% Ut (30% dip in Ut) for 25 cycles <5% Ut (>95% dip in Ut) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Varios 370 / Varios 370 Lux requires continued operation during power mains interruptions, it is recommended that the Varios 370 / Varios 370 Lux be powered from an uninterruptible power supply or a battery.
Power frequency (50/60Hz) magnetic field EN/IEC61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE: Ut is the a.c. mains voltage prior to application of the test level.			

Guidance and manufacturer's declaration - electromagnetic immunity

The Varios 370 / Varios 370 Lux is intended for use in the electromagnetic environment specified below. The customer or the user of the Varios 370 / Varios 370 Lux should assure that it is used in such an environment.

Immunity test	EN/IEC60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF EN/IEC61000-4-6	3Vrms 150 kHz to 80MHz	3Vrms	Portable and mobile RF communications equipment should be used no closer to any part of the Varios 370 / Varios 370 Lux, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 1.2\sqrt{P}$ $d = 1.2\sqrt{P}$ 80MHz to 800MHz $d = 2.3\sqrt{P}$ 800MHz to 2.5GHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF EN/IEC61000-4-3	3V/m 80MHz to 2.5 GHz	3V/m	

NOTE 1 At 80MHz and 800MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobiles radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Varios 370 / Varios 370 Lux is used exceeds the applicable RF compliance level above, the Varios 370 / Varios 370 Lux should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the Varios 370 / Varios 370 Lux.

b Over the frequency range 150kHz to 80MHz, field strengths should be less than 3 V/m.

Cables and accessories	Maximum length	Complies with
Handpiece cord	2 m	RF emissions, CISPR11, EN55011
Foot Control	5 m	Harmonic emissions,
AC Adaptor	5m	Voltage fluctuations/ flicker emission,
		Electrostatic discharge (ESD)
		Electric fast transient / burst
		Surge
		Voltage dips, short interruptions and voltage variations on power supply input lines
		Power frequency(50/60Hz) magnetic field
		Conducted RF
		Radiated RF

Recommended separation distances between portable and mobile RF communications equipment and the Varios 370 / Varios 370 Lux.

The Varios 370 / Varios 370 Lux is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Varios 370 / Varios 370 Lux can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Varios 370 / Varios 370 Lux as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150kHz to 80MHz $d=1.2\sqrt{P}$	80MHz to 800MHz $d=1.2\sqrt{P}$	800MHz to 2,5GHz $d=2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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NSK America Corp www.nskdental.com